



Anti-Allatostatin 1 polyclonal antibody (DPAB-DC3968)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to allatostatin 1. The sequence is CAPSGAQRLYGFGFL
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Allatostatins are neuropeptide hormones in insects and crustacea. They have a twofold function: they both inhibit the generation of juvenile hormone and reduce their food intake. They are therefore putative targets for insecticide research.
Keywords	Allatostatin 1