



Anti-Allopregnanolone polyclonal antibody (DPAB-DC3972)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	5a-PREGNAN-3a-OL-20-ONE HEMISUCCINATE.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Allopregnanolone (3α-hydroxy-5α-pregnan-20-one or 3α,5α-tetrahydroprogesterone), generally abbreviated as ALLO or as 3α,5α-THP, is an endogenous inhibitory pregnane neurosteroid. It is synthesized from progesterone, and is a potent positive allosteric mod
Keywords	ALLO; Allopregnanolone