



Anti-Alpha-Tocopherol polyclonal antibody (DPAB-DC3974)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Alpha-tocopherol acid succinate conjugated with carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ELISA,
Format	Lyophilized
Size	500 µg
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Vitamin E refers to a group of ten lipid-soluble compounds that include both tocopherols and tocotrienols. Of the many different forms of vitamin E, γ-tocopherol is the most common in the North American diet. γ-Tocopherol can be found in corn oil, soybean
---------------------	---

Keywords

Vitamin E; Alpha-Tocopherol (vitamin E)
