



Anti-Amitriptyline polyclonal antibody (DPAB-DC3976)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Amitriptyline conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Amitriptyline (Elavil, Endep, Levate and many others) is a tricyclic antidepressant (TCA). It is the most widely used TCA and is an efficacious treatment for major depressive disorder (clinical depression).
Keywords	Amitriptyline