



Anti-*E. cloacae* ampC Polyclonal antibody (DPAB-DC3977)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>E. cloacae</i> ampC
Immunogen	Native purified <i>Enterobacter cloacae</i> ampC .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>E. cloacae</i>
Conjugate	Unconjugated
Applications	WB, IP, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	AMPC
-----------------	------