



Anti-B. amyloliquefacien amyS Polyclonal antibody (DPAB-DC3980)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	B. amyloliquefacien amyS
Immunogen	Native purified Bacillus amyloliquefaciens amyS.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	B. amyloliquefacien
Conjugate	Unconjugated
Applications	WB, IP, ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot after reconstitution to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Bacillus amyloliquefaciens is a species of bacterium in the genus Bacillus that is the source of the BamH1 restriction enzyme. It also synthesizes a natural antibiotic protein barnase, a widely studied ribonuclease that forms a famously tight complex with
---------------------	---

