



Anti-snake Ancystron (native) polyclonal antibody (DPAB-DC3981)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody has a single antigen reactivity with minimal cross reactivity with other human serum proteins.
Immunogen	Native purified snake (Halys halys) venom Ancystron.
Source/Host	Sheep
Species Reactivity	Human, Snake
Conjugate	Unconjugated
Applications	ICC, IF, ELISA, ID, Immunoelectron Microscopy,
Format	Liquid
Size	1 ml
Buffer	In serum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Ancystron
-----------------	-----------