



Anti-Androstenedione polyclonal antibody (DPAB-DC3983)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Androstenedione(3)-BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Androstenediones are a pair of isomeric steroids which differ at the 5-position. They are distinguished as 5 α -androstenedione and 5 β -androstenedione. In humans, 5 α -androstenedione is a metabolic precursor of both testosterone and estrone.
Keywords	Androstenedione; Androstenedione(3)