



Anti-Anthrax Anthrax Lethal Factor (C-terminal) polyclonal antibody (DPAB-DC3984)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to C-terminus 16 amino acids of Anthrax Lethal Factor.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bacteria
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Anthrax is an acute disease caused by the bacterium <i>Bacillus anthracis</i> . Most forms of the disease are lethal, and it affects both humans and animals. Effective vaccines against anthrax are now available, and some forms of the disease respond well to antib
Keywords	Anthrax Lethal Factor