



Anti-IAV Neuraminidase Polyclonal antibody (DPAB-DC3998)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	IAV Neuraminidase
Immunogen	A synthetic peptide corresponding to C-terminus 15 amino acids of Avian influenza Neuraminidase.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Avian influenza — known informally as avian flu or bird flu — refers to "influenza caused by viruses adapted to birds." The version with the greatest concern is highly pathogenic avian
---------------------	--

influenza (HPAI).

"Bird flu" is a phrase similar to "swine flu," "dog

Keywords

Neuraminidase; Avian Influenza Neuraminidase
