



# Anti-IAV Nonstructural Protein 2 Polyclonal antibody (DPAB-DC4000)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	IAV Nonstructural Protein 2
Immunogen	A synthetic peptide corresponding to internal region 14 amino acids of Avian Influenza NS2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three months. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Avian influenza — known informally as avian flu or bird flu — refers to "influenza caused by viruses adapted to birds." The version with the greatest concern is highly pathogenic avian influenza (HPAI).
--------------	--

"Bird flu" is a phrase similar to "swine flu," "dog

---

**Keywords**

NS2; Avian Influenza NS2

---