



Anti-*L. variegatus* AXIN (C-terminal) polyclonal antibody (DPAB-DC4001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to C-terminus of <i>Lytechinus variegatus</i> AXIN.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Lytechinus variegatus</i>
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<i>Lytechinus variegatus</i> , commonly called the green sea urchin or the variegated sea urchin, is a species of sea urchin that can be found in the warm waters of the western Atlantic Ocean and Caribbean Sea.
Keywords	AXIN