



Anti-Azelaic Acid polyclonal antibody (DPAB-DC4002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Azelaic acid conjugated with carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC-P, IHC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Azelaic acid is an organic compound with the formula (CH ₂) ₇ (CO ₂ H) ₂ . This saturated dicarboxylic acid exists as a white powder. It is found in wheat, rye, and barley. It is a component of a number of hair and skin conditioners.
--------------	--

Keywords

Azelaic acid
