



Anti- Bacitracin Polyclonal antibody (DPAB-DC4004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Bacitracin
Immunogen	Bacitracin conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Bacitracin is a mixture of related cyclic polypeptides produced by organisms of the licheniformis group of <i>Bacillus subtilis</i> var Tracy, first isolated in 1945. These peptides disrupt both gram positive and gram negative bacteria by interfering with cell w
Keywords	Bacitracin