



# Anti-Beta-Defensin-4 (N-terminal) polyclonal antibody (DPAB-DC4010)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus 11 amino acids of Beta-Defensin-4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	ELISA, RIA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from purified antibody in PBS (BSA)
Preservative	None
Storage	Store at -20°C on dry atmosphere. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Beta-Defensin-4
----------	-----------------