



# Anti-Big Gastrin 1 polyclonal antibody (DPAB-DC4015)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Glp-Leu-Gly-Pro-Gln-Gly-Pro-Pro-His-Leu-Val-Ala-Asp-Pro-Ser-Lys-Lys-Gln-Gly-Pro-Trp-Leu-Glu-Glu-Glu-Glu-Glu-Ala-Tyr-Gly-Trp-Met-Asp-Phe-NH <sub>2</sub> conjugate to KLH via carboxyl group.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IP, EIA,
<b>Format</b>	Liquid
<b>Size</b>	150 µg
<b>Buffer</b>	In PBS, pH 7.5 (0.02% sodium azide, 50% glycerol (v/v))
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Big gastrin is a form of gastrin with 34 amino acids in the sequence.
<b>Keywords</b>	Big Gastrin 1