



Anti-M. verrucaria Bilirubin oxidase Polyclonal antibody (DPAB-DC4016)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	M. verrucaria Bilirubin oxidase
Immunogen	Native purified Myrothecium verrucaria Bilirubin oxidase.
Isotype	IgG
Source/Host	Goat
Species Reactivity	M. verrucaria
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IP, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This enzyme belongs to the family of oxidoreductases, to be specific those acting on the CH-CH group of donor with oxygen as acceptor. The systematic name of this enzyme class is
---------------------	--

bilirubin:oxygen oxidoreductase. This enzyme is also called bilirubin oxidase

Keywords

Bilirubin oxidase
