



Anti-Reductase polyclonal antibody (DPAB-DC4017)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to Billverdin Reductase.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Hamster, Human, Mouse, Pig, Rat
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Reduces the gamma-methene bridge of the open tetrapyrrole, biliverdin IX alpha, to bilirubin with the concomitant oxidation of a NADH or NADPH cofactor.
Keywords	Billverdin Reductase