



Anti-S. mansoni 5HT Receptor Polyclonal antibody (DPAB-DC4018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	S. mansoni 5HT Receptor
Immunogen	A synthetic peptide corresponding to C-terminus of Schistosoma mansoni biogenic amine (5HT) receptor.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. mansoni
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The 5-HT receptor (5-hydroxytryptamine receptor) belongs to the G protein coupled receptor 1 family.
---------------------	--

Keywords

5HT Receptor ; biogenic amine (5HT) receptor
