



# Anti-S. mansoni 5HT Receptor Polyclonal antibody (DPAB-DC4019)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                       |
|---------------------------|-------------------------------------------------------------------------------------------------------|
| <b>Target</b>             | S. mansoni 5HT Receptor                                                                               |
| <b>Immunogen</b>          | A synthetic peptide corresponding to N-terminus of Schistosoma mansoni biogenic amine (5HT) receptor. |
| <b>Isotype</b>            | IgG                                                                                                   |
| <b>Source/Host</b>        | Rabbit                                                                                                |
| <b>Species Reactivity</b> | S. mansoni                                                                                            |
| <b>Conjugate</b>          | Unconjugated                                                                                          |
| <b>Applications</b>       | WB, ELISA,                                                                                            |
| <b>Format</b>             | Liquid                                                                                                |
| <b>Size</b>               | 100 µg                                                                                                |
| <b>Buffer</b>             | In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)                                                     |
| <b>Preservative</b>       | 0.01% Sodium Azide                                                                                    |
| <b>Storage</b>            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.   |

## BACKGROUND

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| <b>Introduction</b> | The 5-HT receptor (5-hydroxytryptamine receptor) belongs to the G protein coupled receptor 1 family. |
|---------------------|------------------------------------------------------------------------------------------------------|

**Keywords**

5HT Receptor ; biogenic amine (5HT) receptor

---