



Anti-Boldenone polyclonal antibody (DPAB-DC4021)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Boldenone conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Boldenone (1,4-androstadiene-3-one-17β-ol, available as the undecylenate ester), also known under the trade names Equipoise, Ganabol, Equigan and Ultragan, is an anabolic steroid developed for veterinary use, mostly for treatment of horses. In the US it is
Keywords	Boldenone