



Anti-BRSV Polyclonal antibody (DPAB-DC4024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Bovine Respiratory Syncytial Virus (BRSV) Fusion protein of the RB94 strain.
Target	BRSV
Immunogen	Bovine Respiratory Syncytial Virus (BRSV) Fusion protein of the RB94 strain.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	BRSV
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Bovine Respiratory Syncytial Virus
Keywords	BRSV; Bovine Respiratory Syncytial Virus