



Anti-B. suis VjbR Polyclonal antibody (DPAB-DC4027)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	B. suis VjbR
Immunogen	Recombinant protein corresponding to full length Brucella suis vjbR.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	B. suis
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	VjbR is a LuxR homolog that regulates transcription of many genes.
Keywords	vjbR; Brucella suis vjbR