



Anti-BSA polyclonal antibody (DPAB-DC4028)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA,
Format	Lyophilized
Size	1 mg
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with sterile water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<p>Bovine serum albumin (also known as BSA or "Fraction V") is a serum albumin protein derived from cows. It is often used as a protein concentration standard in lab experiments.</p> <p>The nickname "Fraction V" refers to albumin being the fifth fraction of the ori</p>
Keywords	BSA