



Anti-BSA polyclonal antibody (DPAB-DC4030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Bovine Serum Albumin, BSA.
Source/Host	Sheep
Species Reactivity	Bovine
Conjugate	Unconjugated
Applications	WB (Cell lysate), ICC,
Format	Liquid
Size	200 µg
Buffer	In PBS (0.2% gelatin, 0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store in the dark at 4°C. Do not freeze. Avoid prolonged exposure to light. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Bovine serum albumin (also known as BSA or "Fraction V") is a serum albumin protein derived from cows. It is often used as a protein concentration standard in lab experiments. The nickname "Fraction V" refers to albumin being the fifth fraction of the ori
Keywords	BSA; BSA (FITC)