



# Anti-Buprenorphine polyclonal antibody (DPAB-DC4031)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Buprenorphine conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Buprenorphine is a semi-synthetic, mixed agonist–antagonist opioid receptor modulator that is used to treat opioid addiction in higher dosages, to control moderate acute pain in non-opioid-tolerant individuals in lower dosages and to control moderate chro
Keywords	Buprenorphine