



Anti-B-Zearalanol polyclonal antibody (DPAB-DC4032)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	B-Zearalanol conjugated with BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	1 ml
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	B-Zearalanol
----------	--------------