



Anti-C8 polyclonal antibody (DPAB-DC4033)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-reactivity to other species may exist. Antibody activity to complement components C3-C7, C9 or serum proteins is not present.
Immunogen	Human C8.
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, ELISA, ID,
Format	Liquid
Size	2 ml
Buffer	In 20 mM sodium phosphate buffer, 140 mM NaCl, pH 7.3 (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	C8
-----------------	----