



Anti-Calcitonin polyclonal antibody (DPAB-DC4035)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to human Calcitonin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rabbit, Rat
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Liquid
Size	250 µl
Buffer	In Tris-HCl buffer (1% BSA, 0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Calcitonin (also known as thyrocalcitonin) is a 32-amino acid linear polypeptide hormone that is produced in humans primarily by the parafollicular cells (also known as C-cells) of the thyroid, and in many other animals in the ultimobranchial body. It act
Keywords	Calcitonin

GENE INFORMATION

Entrez Gene ID [796](#)

UniProt ID [P01258](#)
