



# Rabbit anti-cAMP polyclonal antibody (DPAB-DC4038)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |   |
|---------------------|---|
| <b>Immunogen</b>    | A peptide (conjugated to BSA) corresponding to 2-O-succinyl cAMP.                                     |
| <b>Isotype</b>      | IgG   |
| <b>Source/Host</b>  | Rabbit  |
| <b>Conjugate</b>    | Unconjugated  |
| <b>Applications</b> | EIA,  |
| <b>Format</b>       | Liquid  |
| <b>Size</b>         | 100 µl  |
| <b>Buffer</b>       | In PBS (50% glycerol, 1.0% BSA, 0.02% sodium azide)   |
| <b>Preservative</b> | 0.02% Sodium Azide  |
| <b>Storage</b>      | Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|                 |      |
|-----------------|------|
| <b>Keywords</b> | cAMP |
|-----------------|------|