



# Sheep anti-cAMP polyclonal antibody (DPAB-DC4039)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody reacts with cyclic AMP and to lesser degrees with other 5nucleotides.
<b>Immunogen</b>	Succinylated cAMP conjugated with bovine serum albumin.
<b>Source/Host</b>	Sheep
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, RIA,
<b>Format</b>	Liquid
<b>Size</b>	200 µl
<b>Buffer</b>	In serum
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	cAMP
-----------------	------