



# Anti-*C. liberibacter* OMP Polyclonal antibody (DPAB-DC4041)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	<i>C. liberibacter</i> OMP
Immunogen	A synthetic peptide corresponding to internal region of <i>Candidatus Liberibacter asiaticus</i> OMP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>C. liberibacter</i>
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	<i>Candidatus Liberibacter</i> is a genus of gram-negative bacteria in the Rhizobiaceae family. The term <i>Candidatus</i> is used to indicate that it has not proved possible to maintain this bacterium in culture. Detection of the liberibacters is based on PCR amplification.
--------------	---

**Keywords**

C. Liberibacter asiaticus OMP; Candidatus Liberibacter asiaticus OMP

---