



Anti- GluClalpha3A Polyclonal antibody (DPAB-DC4042)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	GluClalpha3A
Immunogen	A synthetic peptide corresponding to C-terminus of <i>Dirofilaria immitis</i> GluClalpha3A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Heartworm (<i>Dirofilaria immitis</i>) is a parasitic roundworm that is spread from host to host through the bites of mosquitoes. The heartworm is a type of filaria, a small thread-like worm, that causes filariasis. The definitive host is the dog, but it can als
--------------	--

Keywords

GluClalpha3A; Canine Heartworm GluClalpha3A
