



Anti- Heartworm Polyclonal antibody (DPAB-DC4043)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Soluble and structural antigens of the adult worm. Antiserum is not absorbed for related filarial Nematodes. Non-reactive with canine serum.
Target	Heartworm
Immunogen	Native purified from <i>Dirofilaria immitis</i> .
Source/Host	Goat
Species Reactivity	Canine
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Heartworm (<i>Dirofilaria immitis</i>) is a parasitic roundworm that is spread from host to host through the bites of mosquitoes. The heartworm is a type of filaria, a small thread-like worm, that
---------------------	--

causes filariasis. The definitive host is the dog, but it can als

Keywords

Canine Heartworm
