



Anti-COPV Envelope Antigen 7 Polyclonal antibody (DPAB-DC4044)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	COPV Envelope Antigen 7
Immunogen	A synthetic peptide corresponding to N-terminus of Canine oral papillomavirus E7 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	COPV
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Papillomaviridae is an ancient taxonomic family of non-enveloped DNA viruses, collectively known as papillomaviruses. Several hundred species of papillomaviruses, traditionally referred to as "types", have been identified infecting all carefully inspected
--------------	---

Keywords

COPV E7; Canine oral papillomavirus E7 protein
