



Anti-COPV Envelope Antigen 6 Polyclonal antibody (DPAB-DC4046)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	COPV Envelope Antigen 6
Immunogen	A synthetic peptide corresponding to C-terminus of Canine oral papillomavirus Protein E6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	COPV
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Papillomaviridae is an ancient taxonomic family of non-enveloped DNA viruses, collectively known as papillomaviruses. Several hundred species of papillomaviruses, traditionally referred to as "types", have been identified infecting all carefully inspected
---------------------	---

Keywords

COPV protein E6; Canine oral papillomavirus Protein E6
