



Anti-Canine papillomavirus type 2 Envelope Antigen 7 Polyclonal antibody (DPAB-DC4047)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Canine papillomavirus type 2 Envelope Antigen 7
Immunogen	A synthetic peptide corresponding to C-terminus of Canine papillomavirus 2 CfpV-2 E7 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Canine papillomavirus
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Papillomaviridae is an ancient taxonomic family of non-enveloped DNA viruses, collectively known as papillomaviruses. Several hundred species of papillomaviruses, traditionally referred to as "types", have been identified infecting all carefully inspected
---------------------	---

Keywords

CfPV-2 E7; Canine papillomavirus type 2 E7 protein
