



# Anti-Canine papillomavirus type 2 Envelope Antigen 7 Polyclonal antibody (DPAB-DC4048)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Canine papillomavirus type 2 Envelope Antigen 7
<b>Immunogen</b>	A synthetic peptide corresponding to N-terminus of Canine papillomavirus 2 CfPV-2 E7 protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Canine papillomavirus
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Papillomaviridae is an ancient taxonomic family of non-enveloped DNA viruses, collectively known as papillomaviruses. Several hundred species of papillomaviruses, traditionally referred to as "types", have been identified infecting all carefully inspected
---------------------	---

**Keywords**

CfPV-2 E7; Canine papillomavirus type 2 E7 protein

---