



Anti-Caproic acid polyclonal antibody (DPAB-DC4051)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Caproic acid conjugated with carrier protein.
Source/Host	Rat
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC-P, IHC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Capric acid occurs naturally in coconut oil (about 10%) and palm kernel oil (about 4%), otherwise it is uncommon in typical seed oils. It is found in the milk of various mammals and to a lesser extent in other animal fats.
Keywords	Capric acid; Caproic acid