



Anti-Caprylic acid polyclonal antibody (DPAB-DC4052)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Caprylic acid conjugated with carrier protein.
Source/Host	Rat
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC-P, IHC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Caprylic acid is the common name for the eight-carbon saturated fatty acid known by the systematic name octanoic acid. Its compounds are found naturally in the milk of various mammals, and as a minor constituent of coconut oil and palm kernel oil. It is a
Keywords	Caprylic acid