



Anti-Casomorphin-beta polyclonal antibody (DPAB-DC4057)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to human Casomorphin-beta. The sequence is CYPFVEPI
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Casomorphins are peptides, i.e., protein fragments, derived from the digestion of milk protein casein. The distinguishing characteristic of casomorphins is that they have an opioid effect. β-casomorphins are found in cheeses made from bovine milk; their c
Keywords	Casomorphin-beta