



Anti- CASP1 Polyclonal antibody (DPAB-DC4058)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CASP1
Immunogen	A synthetic peptide corresponding to N-terminus of fall armyworm CASP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fall armyworm
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CASP1
----------	-------

GENE INFORMATION

Entrez Gene ID [834](#)

UniProt ID [P29466](#)
