



Anti-CASP13 (aa 91-105) polyclonal antibody (DPAB-DC4059)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 91-105 of human CASP13.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CASP13
----------	--------

GENE INFORMATION

Entrez Gene ID [837](#)

UniProt ID [P49662](#)
