



Anti-Cav1b (N-terminal) polyclonal antibody (DPAB-DC4064)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of <i>Lymnaea stagnalis</i> Cav1b.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Sea urchin
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<i>Lymnaea stagnalis</i> , better known as the great pond snail, is a species of large air-breathing freshwater snail, an aquatic pulmonate gastropod mollusk in the family Lymnaeidae.
Keywords	Cav1b