



Anti-CCL3L1 polyclonal antibody (DPAB-DC4068)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to 70 amino acids of swine CCL3L1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Pig
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Chemokine (C-C motif) ligand 3-like 1, also known as CCL3L1, is a protein which in humans is encoded by the CCL3L1 gene
Keywords	CCL3L1

GENE INFORMATION

Entrez Gene ID [6349](#)

UniProt ID [P16619](#)
