



# Anti-CCR5 polyclonal antibody (DPAB-DC4069)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody is specific to CCR5. It demonstrates 2D7-like reactivity.
<b>Immunogen</b>	A synthetic peptide corresponding to CCR5. The sequence is EWQKEGLVTLWLGGGGSC
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Primates
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC, ELISA, Flow Cyt,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In antiserum
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	CCR5
-----------------	------

## GENE INFORMATION

Entrez Gene ID [1234](#)

---

UniProt ID [P51681](#)

---