



# Anti-Cell Hemolysin (native) polyclonal antibody (DPAB-DC4075)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Native from sheep washed pooled red blood cells (RBC).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Sheep
Conjugate	Unconjugated
Applications	Agg,
Format	Lyophilized
Size	2 ml
Buffer	Lyophilized from antiserum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Cell Hemolysin
----------	----------------