



Rabbit anti-cGMP polyclonal antibody (DPAB-DC4076)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A peptide corresponding to 2O-succinyl guanosine-3,5 cyclic monophosphate (cGMP).
Isotype	IgG
Source/Host	Rabbit
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In PBS (50% glycerol, 1.0% BSA, 0.02% thimerosal)
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	cGMP
-----------------	------