



Anti-Chlamydia LGV2 Polyclonal antibody (DPAB-DC4077)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts predominantly against LPS antigens with some activity against MOMP by Western blot. Cross-reacts against all <i>C. trachomatis</i> serovars, <i>C. psittaci</i> and <i>C. pneumoniae</i> .
Target	Chlamydia LGV2
Immunogen	Elementary bodies from LGV2 serovar.
Source/Host	Pig
Species Reactivity	Chlamydia
Conjugate	Unconjugated
Applications	IF, ELISA, Capture Antibody, Detection Antibody,
Format	Liquid
Size	1 ml
Buffer	In antiserum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords Chlamydia; Chlamydia (Genus Specific)