



Anti-Chlamydia LGV2 Polyclonal antibody (DPAB-DC4078)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Predominant reactivity against LPS antigen with some MOMP by Western blot. Cross-reacts against all <i>C. trachomatis</i> serovars, <i>C. psittaci</i> and <i>C. pneumoniae</i> .
Target	Chlamydia LGV2
Immunogen	Elementary bodies from LGV2 serovar.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chlamydia
Conjugate	Unconjugated
Applications	WB, ELISA, Detection Antibody,
Format	Liquid
Size	1 ml
Buffer	In buffer containing proclin 300
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords

Chlamydia; Chlamydia (Genus Specific)
